

Submitting clinician:				
Practice name & address:				
Tel:	Fax:			
E-mail:				
Owner name & address:	Animal Name:	I.D.:	Breed:	
	Sex:	Age:		
Preferred reporting method: Fax E-mail Date collected:				
 Sampling instructions: The test is performed on the 1st and 3rd Thursdays of each month only. A minimum of 3ml of heparinised blood should be collected on the Tuesday or Wednesday before the test is run and sent to Biobest by an overnight service such as Royal Mail Special Delivery Guaranteed. Samples must not be chilled as this will interfere with cell viability and test performance. 				
Test Interpretation and Additional Guidance: This test is intended to provide support for a diagnosis of canine mycobacterial disease and whether the infection is likely to be due to bovine TB (bTB) or other Mycobacterial species. If suspected at <i>PM</i> bTB is a notifiable disease The test takes a minimum 96 hours to complete therefore results can be expected on the Tuesday following set up. Biobest may report summary test results to the TB team at APHA.				
Clinical History:				

Biobest Use Only			
Date of Receipt:	Databased:	Request Form Number:	
Reported:	QC:	Report Faxed/Emailed:	

Biobest Laboratories Ltd, 6 Charles Darwin House, The Edinburgh Technopole, Milton Bridge, Nr Penicuik, EH26 0PY, UK Tel: +44 (0)131 440 2628 Fax: +44 (0)131 440 9587 email: enquiry@biobest.co.uk www.biobest.co.uk