

## CANINE SUBMISSION FORM Interferon gamma (Y) testing for suspected tuberculosis cases

Submitting clinician:					
Practice name & address:					
Tel:	Fax	Fax:			
E-mail:					
Owner name & address:	Ani Nar	mal ne:	I.D.:	Breed:	
	Sex	:	Age:		
Preferred reporting method: Fax   E-mail   Date collected:			:		
Sampling instructions: The test is performed every Thursday. A minimum of 3ml of heparinised blood should be collected on the Tuesday before the test is run and sent to Biobest by an overnight service such as Royal Mail Special Delivery Guaranteed.  Samples must not be chilled as this will interfere with cell viability and test performance.  Test Interpretation and Additional Guidance:					
This test is intended to provide support for a diagnosis of canine mycobacterial disease and whether the infection is likely to be due to bovine TB (bTB) or other Mycobacterial species. If suspected at <i>PM</i> bTB is a notifiable disease The test takes a minimum 96 hours to complete therefore results can be expected on the Tuesday following set up. Biobest may report summary test results to the TB team at APHA.					
Clinical History:					
Disharat Hara Onda					
Date of Receipt:	Biobest Use Only  Databased:		Paguest Form Nu	Request Form Number:	
			·		
Reported:	QC:	Report Faxed/Emailed:			

Biobest Laboratories Ltd, 6 Charles Darwin House, The Edinburgh Technopole, Milton Bridge, Nr Penicuik, EH26 0PY, UK Tel: +44 (0)131 440 2628 Fax: +44 (0)131 440 9587 email: enquiry@biobest.co.uk www.biobest.co.uk